

兔抗 RCOR3 多克隆抗体

中文名称: 兔抗 RCOR3 多克隆抗体

英文名称: Anti-RCOR3 rabbit polyclonal antibody

别名: REST corepressor 3

抗原: RCOR3

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	May act as a component of a corepressor complex that represses transcription.
Applications:	ELISA, IHC
Name of antibody:	RCOR3
Immunogen:	Fusion protein of human RCOR3
Full name:	REST corepressor 3
SwissProt:	Q9P2K3
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human lung cancer and Human esophagus cancer
IHC Recommend dilution:	50-300

