

兔抗 RILPL2 多克隆抗体

中文名称：兔抗 RILPL2 多克隆抗体

英文名称：Anti-RILPL2 rabbit polyclonal antibody

别名：RLP2

相关类别：一抗

抗原：RILPL2

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a protein that contains a rab-interacting lysosomal protein-like domain. This protein may be involved in regulating lysosome morphology. This protein may also be a target for the Hepatitis C virus and assist in viral replication. Alternate splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	RILPL2
Immunogen:	Full length fusion protein
Full name:	Rab interacting lysosomal protein-like 2
Synonyms :	RLP2

SwissProt:	Q969X0
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human brain and Human prostate cancer
IHC Recommend dilution:	25-100

