

兔抗 RMND5A 多克隆抗体

- 中文名称: 兔抗 RMND5A 多克隆抗体
- 英文名称: Anti-RMND5A rabbit polyclonal antibody
- 别名: required for meiotic nuclear division 5 homolog A; CTLH; GID2; RMD5; GID2A; p44CTLH
- 储存: 冷冻 (-20℃)
- 抗原: RMND5A
- 宿主: Rabbit
- 反应种属: Human, Mouse
- 相关类别: 一抗
- 标记物: Unconjugate
- 克隆类型: rabbit polyclonal

技术规格

Background:	The function of this protein remains unknown.
Applications:	ELISA, IHC
Name of antibody:	RMND5A
Immunogen:	Fusion protein of human RMND5A
Full name:	required for meiotic nuclear division 5 homolog A
Synonyms:	CTLH; GID2; RMD5; GID2A; p44CTLH
SwissProt:	Q9H871
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer
IHC Recommend dilution:	20-100

