

## 兔抗 RNF114 多克隆抗体

中文名称：兔抗 RNF114 多克隆抗体

英文名称：Anti-RNF114 rabbit polyclonal antibody

别名：ZNF313; PSORS12

抗原：RNF114

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	E3 ubiquitin-protein ligase promoting the ubiquitination and degradation of the CDK inhibitor CDKN1A and probably also CDKN1B and CDKN1C. These activities stimulate cell cycle's G1-to-S phase transition and suppress cellular senescence. May play a role in spermatogenesis.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	RNF114
<b>Immunogen:</b>	Full length fusion protein
<b>Full name:</b>	ring finger protein 114
<b>Synonyms :</b>	ZNF313; PSORS12
<b>SwissProt:</b>	Q9Y508

<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human prostate cancer and Human esophagus cancer
<b>IHC Recommend dilution:</b>	50-100

