

## 兔抗 RNF115 多克隆抗体

中文名称：兔抗 RNF115 多克隆抗体

英文名称：Anti-RNF115 rabbit polyclonal antibody

别名：BCA2; ZNF364

储存：冷冻 (-20℃)

抗原：RNF115

宿主：Rabbit

相关类别：一抗

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	E3 ubiquitin-protein ligase that mediates E2-dependent, 'Lys-48'- and/or 'Lys-63'-linked polyubiquitination of substrates and may play a role in diverse biological processes. Through their polyubiquitination, may play a role in the endosomal trafficking and degradation of membrane receptors including EGFR, FLT3, MET and CXCR4.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	RNF115
<b>Immunogen:</b>	Full length fusion protein
<b>Full name:</b>	ring finger protein 115
<b>Synonyms :</b>	BCA2; ZNF364

<b>SwissProt:</b>	Q9Y4L5
<b>ELISA Recommended dilution:</b>	2000-10000
<b>IHC positive control:</b>	Human liver cancer and Human prostate cancer
<b>IHC Recommend dilution:</b>	30-150

