

兔抗 RPF2 多克隆抗体

中文名称：兔抗 RPF2 多克隆抗体

英文名称：Anti-RPF2 rabbit polyclonal antibody

别名：BXDC1; bA397G5.4

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：RPF2

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

BXDC1 (Brix domain-containing protein 1) is a 306 amino acid protein encoded by the human gene BXDC1. BXDC1 is a nuclear protein that contains one Brix domain. Brix domain containing proteins represent a family of proteins involved in the biogenesis of large ribosomal subunit

	<p>s. The Brix domain is a region with homology to the yeast protein Pitx1 (Ribosome biogenesis protein BRX1). Pitx1 is part of a complex that includes RPF1, RPF2 and SSF1 or SSF2. This complex is required for the biogenesis of the 60S ribosomal subunit.</p>
Applications:	ELISA, WB
Name of antibody:	RPF2
Immunogen:	Synthetic peptide of human RPF2
Full name:	ribosome production factor 2 homolog
Synonyms:	BXDC1; bA397G5.4
SwissProt:	Q9H7B2
ELISA Recommended dilution:	5000-10000
WB Predicted band size:	36 kDa
WB Positive control:	K562,Hela and A549 cell lysates
WB Recommended dilution:	200-1000

