

兔抗 RRP1 多克隆抗体

中文名称：兔抗 RRP1 多克隆抗体

英文名称：Anti-RRP1 rabbit polyclonal antibody

别名：ribosomal RNA processing 1; NNP-1; NOP52; RRP1A; D21S2056E

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：RRP1

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is the putative homolog of the yeast ribosomal RNA processing protein RRP1. The encoded protein is involved in the late stages of nucleogenesis at the end of mitosis, and may be required for the generation of 28S rRNA.
Applications:	ELISA, WB, IHC
Name of antibody:	RRP1
Immunogen:	Fusion protein of human RRP1
Full name:	ribosomal RNA processing 1
Synonyms:	NNP-1; NOP52; RRP1A; D21S2056E
SwissProt:	P56182
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human lung cancer
IHC Recommend dilution:	40-200
WB Predicted band size:	53 kDa
WB Positive control:	HepG2, HeLa, 293T and Jurkat cell lysates
WB Recommended dilution:	500-2000



