

兔抗 SAXO2 多克隆抗体

中文名称：兔抗 SAXO2 多克隆抗体

英文名称：Anti-SAXO2 rabbit polyclonal antibody

别名：stabilizer of axonemal microtubules 2; FAM154B

抗原：SAXO2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The function of this protein remains unknown.
Applications:	ELISA, IHC
Name of antibody:	SAXO2
Immunogen:	Synthetic peptide of human SAXO2
Full name:	stabilizer of axonemal microtubules 2
Synonyms:	FAM154B
SwissProt:	Q658L1
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	40-200

