

**兔抗 SEC11A 多克隆抗体**

中文名称：兔抗 SEC11A 多克隆抗体

英文名称：Anti-SEC11A rabbit polyclonal antibody

别名：SEC11 homolog A, signal peptidase complex subunit; SPC18; SPCS4A; SEC11L1; sid2895; 1810012E07Rik

相关类别：一抗

储存：冷冻 (-20℃)

抗原：SEC11A

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

**技术规格**

<b>Background:</b>	This gene encodes a member of the peptidase S26B family. The encoded protein is an 18kDa subunit of the signal peptidase complex and has been linked to cell migration and invasion, gastric cancer and lymph node metastasis. Alternative splicing results in multiple transcript variants. A related pseudogene has been identified on chromosome 8.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	SEC11A
<b>Immunogen:</b>	Fusion protein of human SEC11A
<b>Full name:</b>	SEC11 homolog A, signal peptidase complex subu

	nit
<b>Synonyms:</b>	SPC18; SPCS4A; SEC11L1; sid2895; 1810012E07Rik
<b>SwissProt:</b>	P67812
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human gastric cancer
<b>IHC Recommend dilution:</b>	40-200
<b>WB Predicted band size:</b>	21 kDa
<b>WB Positive control:</b>	TM4, HEPG2 and A549 cell, Mouse liver tissue lysates
<b>WB Recommended dilution:</b>	500-2000



