

## 兔抗 SEC22B 多克隆抗体

中文名称：兔抗 SEC22B 多克隆抗体

英文名称： Anti-SEC22B rabbit polyclonal antibody

别名： SEC22 homolog B, vesicle trafficking protein (gene/pseudogene); ERS-24; SEC22L1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SEC22B

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a member of the SEC22 family of vesicle trafficking proteins. It seems to complex with SNARE and it is thought to play a role in the ER-Golgi protein trafficking. This protein has strong similarity to Mus musculus and Cricetulus griseus proteins.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	SEC22B
<b>Immunogen:</b>	Fusion protein of human SEC22B
<b>Full name:</b>	SEC22 homolog B, vesicle trafficking protein (gene/pseudogene)
<b>Synonyms:</b>	ERS-24; SEC22L1

<b>SwissProt:</b>	O75396
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human tonsil
<b>IHC Recommend dilution:</b>	50-300
<b>WB Predicted band size:</b>	25 kDa
<b>WB Positive control:</b>	A172, LOVO and Hela cell lysates
<b>WB Recommended dilution:</b>	500-2000



