

兔抗 SERPINA1 多克隆抗体

中文名称：兔抗 SERPINA1 多克隆抗体

英文名称：Anti-SERPINA1 rabbit polyclonal antibody

别名：PI; A1A; AAT; PI1; A1AT; PRO2275; alpha1AT

相关类别：一抗

抗原：SERPINA1

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

| | |
|--------------------------|--|
| Background: | The protein encoded by this gene is secreted and is a serine protease inhibitor whose targets include elastase, plasmin, thrombin, trypsin, chymotrypsin, and plasminogen activator. Defects in this gene can cause emphysema or liver disease. Several transcript variants encoding the same protein have been found for this gene. |
| Applications: | ELISA, IHC |
| Name of antibody: | SERPINA1 |
| Immunogen: | Fusion protein of human SERPINA1 |
| Full name: | serpin peptidase inhibitor, clade A (alpha-1 antiprotease, antitrypsin), member 1 |
| Synonyms : | PI; A1A; AAT; PI1; A1AT; PRO2275; alpha1AT |

| | |
|------------------------------------|-------------|
| SwissProt: | P01009 |
| ELISA Recommended dilution: | 2000-5000 |
| IHC positive control: | Human brain |
| IHC Recommend dilution: | 50-200 |

