

兔抗 SERTAD1 多克隆抗体

中文名称：兔抗 SERTAD1 多克隆抗体

英文名称：Anti-SERTAD1 rabbit polyclonal antibody

别名：SEI1; TRIP-Br1

储存：冷冻（-20℃）

抗原：SERTAD1

宿主：Rabbit

相关类别：一抗

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

background:	Acts at E2F-responsive promoters as coregulator to integrate signals provided by PHD- and/or bromodomain-containing transcription factors. Stimulates E2F1/TFDP1 transcriptional activity. Renders the activity of cyclin D1/CDK4 resistant to the inhibitory effects of CDKN2A/p16INK4A.
Applications:	ELISA, IHC
Name of antibody:	SERTAD1
Immunogen:	Full length fusion protein
Full name:	SERTA domain containing 1
Synonyms:	SEI1; TRIP-Br1
SwissProt:	Q9UHV2

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer and human cervical cancer
IHC Recommend dilution:	25-100

