

兔抗 SESN3 多克隆抗体

中文名称: 兔抗 SESN3 多克隆抗体

英文名称: Anti-SESN3 rabbit polyclonal antibody

别名: SEST3

储存: 冷冻 (-20℃)

抗原: SESN3

宿主: Rabbit

反应种属: Human, Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a member of the sestrin family of stress-induced proteins. The encoded protein reduces the levels of intracellular reactive oxygen species induced by activated Ras downstream of RAC-alpha serine/threonine-protein kinase (Akt) and FoxO transcription factor. The protein is required for normal regulation of blood glucose, insulin resistance and plays a role in lipid storage in obesity. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, WB
Name of antibody:	SESN3
Immunogen:	Full length fusion protein
Full name:	sestrin 3

Synonyms :	SEST3
SwissProt:	P58005
ELISA Recommended dilution:	1000-2000
WB Predicted band size:	57 kDa
WB Positive control:	Mouse heart tissue
WB Recommended dilution:	200-1000

