

## 兔抗 SFN 多克隆抗体

中文名称：兔抗 SFN 多克隆抗体

英文名称： Anti-SFN rabbit polyclonal antibody

别名： YWHAS

抗原： SFN

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	Adapter protein implicated in the regulation of a large spectrum of both general and specialized signaling pathways. Binds to a large number of partners, usually by recognition of a phosphoserine or phosphothreonine motif. Binding generally results
--------------------	--

	<p>in the modulation of the activity of the binding partner. When bound to KRT17, regulates protein synthesis and epithelial cell growth by stimulating Akt/mTOR pathway. May also regulate MDM2 autoubiquitination and degradation and thereby activate p53/TP53.</p>
<b>Applications:</b>	ELISA,WB, IHC
<b>Name of antibody:</b>	SFN
<b>Immunogen:</b>	Full length fusion protein
<b>Full name:</b>	stratifin
<b>Synonyms :</b>	YWHAS
<b>SwissProt:</b>	P31947
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer and Human breast cancer
<b>IHC Recommend dilution:</b>	500-2000

