

## 兔抗 SGCD 多克隆抗体

中文名称：兔抗 SGCD 多克隆抗体

英文名称： Anti-SGCD rabbit polyclonal antibody

别名： sarcoglycan delta; SGD; DAGD; 35DAG; CMD1L; SGCDP; SG-delta

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SGCD

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is one of the four known components of the sarcoglycan complex, which is a subcomplex of the dystrophin-glycoprotein complex (DGC). DGC forms a link between the F-actin cytoskeleton and the extracellular matrix. This protein is expressed most abundantly in skeletal and cardiac muscle. Mutations in this gene have been associated with autosomal recessive limb-girdle muscular dystrophy and dilated cardiomyopathy. Alternatively spliced transcript variants encoding distinct isoforms have been observed for this gene.
<b>Applications:</b>	ELISA, WB, IHC

<b>Name of antibody:</b>	SGCD
<b>Immunogen:</b>	Fusion protein of human SGCD
<b>Full name:</b>	sarcoglycan delta
<b>Synonyms:</b>	SGD; DAGD; 35DAG; CMD1L; SGCDP; SG-delta
<b>SwissProt:</b>	Q92629
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human tonsil
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	32 kDa
<b>WB Positive control:</b>	Jurkat cell lysate
<b>WB Recommended dilution:</b>	200-1000



