

## 兔抗 SGK2 多克隆抗体

中文名称: 兔抗 SGK2 多克隆抗体

英文名称: Anti-SGK2 rabbit polyclonal antibody

别名: H-SGK2; dJ138B7.2

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: SGK2

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a serine/threonine protein kinase. Although this gene product is similar to serum- and glucocorticoid-induced protein kinase (SGK), this gene is not induced by serum or glucocorticoids. This gene is induced in response to signals that activate phosphatidylinositol 3-kinase, which is also true for SGK. Alternative splicing results in multiple transcript variants
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	SGK2
<b>Immunogen:</b>	Synthetic peptide of human SGK2
<b>Full name:</b>	serum/glucocorticoid regulated kinase 2

<b>Synonyms:</b>	H-SGK2; dJ138B7.2
<b>SwissProt:</b>	Q9HBY8
<b>IHC positive control:</b>	Human gastric cancer and human thyroid cancer
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	48 kDa
<b>WB Positive control:</b>	Mouse liver tissue
<b>WB Recommended dilution:</b>	200-1000



