

兔抗 SH3GLB2 多克隆抗体

中文名称：兔抗 SH3GLB2 多克隆抗体

英文名称： Anti-SH3GLB2 rabbit polyclonal antibody

别名： RRIG1; PP6569; PP9455

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： SH3GLB2

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

| | |
|--------------------------------|---|
| Background: | Belongs to the endophilin family. 2 isoforms of the human protein are produced by alternative splicing. |
| Applications: | WB, IHC |
| Name of antibody: | SH3GLB2 |
| Immunogen: | Synthesized peptide derived from internal of human SH3GLB2. |
| Full name: | SH3-domain GRB2-like endophilin B2 |
| Synonyms : | RRIG1; PP6569; PP9455 |
| SwissProt: | Q9NR46 |
| IHC positive control: | Human prostate carcinoma tissue |
| IHC Recommend dilution: | 50-100 |

| | |
|---------------------------------|------------------|
| WB Predicted band size: | 44 kDa |
| WB Positive control: | COS cells lysate |
| WB Recommended dilution: | 500-3000 |

