

## 兔抗 SIN3B 多克隆抗体

中文名称：兔抗 SIN3B 多克隆抗体

英文名称：Anti-SIN3B rabbit polyclonal antibody

抗原：SIN3B

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

相关类别：一抗

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Acts as a transcriptional repressor. Interacts with MXI1 to repress MYC responsive genes and antagonize MYC oncogenic activities. Interacts with MAD-MAX heterodimers by binding to MAD. The heterodimer then represses transcription by tethering SIN3B to DNA. Also forms a complex with FOXK1 which represses transcription (By similarity).
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	SIN3B
<b>Immunogen:</b>	Synthetic peptide of human SIN3B
<b>Full name:</b>	SIN3 transcription regulator family member B
<b>SwissProt:</b>	O75182

<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human thyroid cancer
<b>IHC Recommend dilution:</b>	25-100

