

## 兔抗 MEF2C 多克隆抗体

中文名称：兔抗 MEF2C 多克隆抗体

英文名称：Anti-MEF2C rabbit polyclonal antibody

别名：DEL5q14.3; C5DELq14.3

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：MEF2C

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This locus encodes a member of the MADS box transcription enhancer factor 2 (MEF2) family of proteins, which play a role in myogenesis. The encoded protein, MEF2 polypeptide C, has both trans-activating and DNA binding activities. This protein may play a role in maintaining the differentiated state of muscle cells. Mutations and deletions at this locus have been associated with severe mental retardation, stereotypic movements, epilepsy, and cerebral malformation. Alternatively spliced transcript variants have been described.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	MEF2C
<b>Immunogen:</b>	Synthetic peptide of human MEF2C
<b>Full name:</b>	myocyte enhancer factor 2C

<b>Synonyms:</b>	DEL5q14.3; C5DELq14.3
<b>SwissProt:</b>	Q06413
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human lung cancer
<b>IHC Recommend dilution:</b>	20-100



