

兔抗 MOB1B 多克隆抗体

中文名称：兔抗 MOB1B 多克隆抗体

英文名称：Anti-MOBKL1A rabbit polyclonal antibody

别名：MATS2; MOB4A; MOBKL1A

相关类别：一抗

抗原：MOB1B

储存：冷冻 (-20℃)

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is similar to the yeast Mob1 protein. Yeast Mob1 binds Mps1p, a protein kinase essential for spindle pole body duplication and mitotic checkpoint regulation. Three transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	MOB1B
Immunogen:	Full length fusion protein
Full name:	MOB kinase activator 1B
Synonyms :	MATS2; MOB4A; MOBKL1A
SwissProt:	Q7L9L4
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human thyroid cancer and Human ovarian cancer
IHC Recommend dilution:	25-100



