

兔抗 MORF4L1 多克隆抗体

中文名称：兔抗 MORF4L1 多克隆抗体

英文名称：Anti-MORF4L1 rabbit polyclonal antibody

别名：Eaf3; MEAF3; MRG15; FWP006; S863-6; HsT17725; MORFRG15

相关类别：一抗

抗原：MORF4L1

储存：冷冻 (-20℃)

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

The members of the mortality factor family include mortality factor 4 (MORF4), MORF4L1 (also known as MRG15) and MORF4-related gene X (MRGX). The human MORF4 gene maps to chromosome 4q33-q34.1. MORF4 induces a senescent-like phenotype in complementation group B immortal cell lines. The genes encoding MRG15 and MRGX map to chromosomes 15q24 and Xq22, respectively. MORF4, MORF4L1 and MRGX each contain a C-terminal leucine zipper. An association between MORF4L1, Rb (retinoblastoma tumor suppressor) and PAM14 (protein associated with MORF4L1) suggests a role for MORF4L1 in transcription regulation. MORF4L1 also associates with the histo

	ne acetyl transferase MOF. In addition, MORF4 and MRGX interact with mSin3A and TLE (transducin-like enhancer of split). The MORF/ mSin3A/TLE association may repress transcription. In Purkinje cells, MORF4L1 localizes to the dendrites and the nuclei.
Applications:	ELISA, IHC
Name of antibody:	MORF4L1
Immunogen:	Synthetic peptide of human MORF4L1
Full name:	mortality factor 4 like 1
Synonyms :	Eaf3; MEAF3; MRG15; FWP006; S863-6; HsT17725; MORFRG15
SwissProt:	Q9UBU8
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human lung cancer and Human tonsil
IHC Recommend dilution:	30-150

