

## 兔抗 MST1R 多克隆抗体

中文名称：兔抗 MST1R 多克隆抗体

英文名称：Anti-MST1R rabbit polyclonal antibody

别名：macrophage stimulating 1 receptor; RON; PTK8; CD136; CDw136

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：MST1R

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

**Background:**

This gene encodes a cell surface receptor for macrophage-stimulating protein (MSP) with tyrosine kinase activity. The mature form of this protein is a heterodimer of disulfide-linked alpha and beta subunits, generated by proteolytic cleavage of a single-chain precursor. The beta subunit undergoes tyrosine phosphorylation upon stimulation by MSP. This protein is expressed on the ciliated epithelia of the mucociliary transport apparatus of the lung, and together with MSP, thought to be involved in host defense. Alternative splicing generates multiple transcript variants encoding different isoforms that may undergo similar proteolytic processing.

<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	MST1R
<b>Immunogen:</b>	Synthetic peptide of human MST1R
<b>Full name:</b>	macrophage stimulating 1 receptor
<b>Synonyms:</b>	RON; PTK8; CD136; CDw136
<b>SwissProt:</b>	Q04912
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	25-100

