

兔抗 MTHFS 多克隆抗体

中文名称: 兔抗 MTHFS 多克隆抗体

英文名称: Anti-MTHFS rabbit polyclonal antibody

别名: methenyltetrahydrofolate synthetase; HsT19268

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: MTHFS

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is an enzyme that catalyzes the conversion of 5-formyltetrahydrofolate to 5,10-methenyltetrahydrofolate, a precursor of reduced folates involved in 1-carbon metabolism. An increased activity of the encoded protein can result in an increased folate turnover rate and folate depletion. Three transcript variants encoding two different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	MTHFS
Immunogen:	Fusion protein of human MTHFS
Full name:	methenyltetrahydrofolate synthetase
Synonyms:	HsT19268

SwissProt:	P49914
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human prostate cancer and Human brain
IHC Recommend dilution:	40-200



