

兔抗 MUC13 多克隆抗体

中文名称: 兔抗 MUC13 多克隆抗体

英文名称: Anti-MUC13 rabbit polyclonal antibody

别名: DRCC1, MUC-13

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: MUC13

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Synthetic peptide of human MUC13
Applications:	ELISA, IHC
Name of antibody:	MUC13
Immunogen:	Synthetic peptide of human MUC13
Full name:	mucin 13, cell surface associated
Synonyms:	DRCC1, MUC-13
SwissProt:	Q9H3R2
ELISA Recommended dilution:	1000-2000
Positive control:	Human liver cancer;
Recommended dilution:	15-50

