

## 兔抗 SLC24A5 多克隆抗体

中文名称：兔抗 SLC24A5 多克隆抗体

英文名称：Anti-SLC24A5 rabbit polyclonal antibody

级别：JSX; OCA6; NCKX5; SHEP4

抗原：SLC24A5

储存：冷冻（-20℃）

宿主：Rabbit

相关类别：一抗

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene is a member of the potassium-dependent sodium/calcium exchanger family and encodes an intracellular membrane protein with 2 large hydrophilic loops and 2 sets of multiple transmembrane-spanning segments. Sequence variation in this gene has been associated with differences in skin pigmentation.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	SLC24A5
<b>Immunogen:</b>	Synthetic peptide of human SLC24A5
<b>Full name:</b>	solute carrier family 24 (sodium/potassium/calcium exchanger), member 5
<b>Synonyms:</b>	JSX; OCA6; NCKX5; SHEP4

<b>SwissProt:</b>	Q71RS6
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human thyroid cancer
<b>IHC Recommend dilution:</b>	30-150

