

兔抗 SLC38A1 多克隆抗体

中文名称：兔抗 SLC38A1 多克隆抗体

英文名称：Anti-SLC38A1 rabbit polyclonal antibody

级别：ATA1; NAT2; SAT1; SNAT1

抗原：SLC38A1

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

相关类别：一抗

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Amino acid transporters play essential roles in the uptake of nutrients, production of energy, chemical metabolism, detoxification, and neurotransmitter cycling. SLC38A1 is an important transporter of glutamine, an intermediate in the detoxification of ammonia and the production of urea. Glutamine serves as a precursor for the synaptic transmitter, glutamate.
Applications:	ELISA, WB
Name of antibody:	SLC38A1
Immunogen:	Synthetic peptide of human SLC38A1
Full name:	solute carrier family 38 member 1

Synonyms:	ATA1; NAT2; SAT1; SNAT1
SwissProt:	Q9H2H9
ELISA Recommended dilution:	5000-10000
WB Predicted band size:	54 kDa
WB Positive control:	Human liver cancer tissue lysate
WB Recommended dilution:	200-1000

