

兔抗 SLC8A3 多克隆抗体

中文名称：兔抗 SLC8A3 多克隆抗体

英文名称：Anti-SLC8A3 rabbit polyclonal antibody

别名：NCX3

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：SLC8A3

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the sodium/calcium exchanger integral membrane protein family. Three mammalian isoforms in family 8 have been identified. Na ⁺ /Ca ²⁺ exchange proteins are involved in maintaining Ca ²⁺ homeostasis in a wide variety of cell types. The protein is regulated by intracellular calcium ions and is found in both the plasma membrane and intracellular organellar membranes, where exchange of Na ⁺ for Ca ²⁺ occurs in an electrogenic manner. Alternative splicing has been observed for this gene and multiple variants have been described.
Applications:	ELISA, IHC
Name of antibody:	SLC8A3

Immunogen:	Synthetic peptide of human SLC8A3
Full name:	Solute carrier family 8 (sodium/calcium exchanger), member 3
Synonyms :	NCX3
SwissProt:	P57103
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human ovarian cancer and Human lung cancer
IHC Recommend dilution:	25-100

