

兔抗 NARS1 多克隆抗体

中文名称: 兔抗 NARS1 多克隆抗体

英文名称: Anti-NARS1 rabbit polyclonal antibody

别名: asparaginyl-tRNA synthetase 1; NARS; ASNRS

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: NARS1

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Aminoacyl-tRNA synthetases are a class of enzymes that charge tRNAs with their cognate amino acids. Asparaginyl-tRNA synthetase is localized to the cytoplasm and belongs to the class II family of tRNA synthetases. The N-terminal domain represents the signature sequence for the eukaryotic asparaginyl-tRNA synthetases. [provided by RefSeq, Jul 2008]
Applications:	ELISA, IHC
Name of antibody:	NARS1
Immunogen:	Fusion protein of human NARS1
Full name:	asparaginyl-tRNA synthetase 1
Synonyms:	NARS; ASNRS
SwissProt:	O43776

IHC positive control:	Human tonsil
IHC Recommend dilution:	40-200



