

## 兔抗 NCBP1 多克隆抗体

中文名称: 兔抗 NCBP1 多克隆抗体

英文名称: Anti-NCBP1 rabbit polyclonal antibody

别名: nuclear cap binding protein subunit 1; NCBP; Sto1; CBP80

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: NCBP1

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	The product of this gene is a component of the nuclear cap-binding protein complex (CBC), which binds to the monomethylated 5' cap of nascent pre-mRNA in the nucleoplasm. The encoded protein promotes high-affinity mRNA -cap binding and associates with the CTD of RNA polymerase II. The CBC promotes pre-mRNA splicing, 3'-end processing, RNA nuclear export, and nonsense-mediated mRNA decay.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	NCBP1
<b>Immunogen:</b>	Synthetic peptide of human NCBP1
<b>Full name:</b>	nuclear cap binding protein subunit 1

<b>Synonyms:</b>	NCBP; Sto1; CBP80
<b>SwissProt:</b>	Q09161
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human lung cancer and Human liver cancer
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	92 kDa
<b>WB Positive control:</b>	Jurkat , NIH/3T3 , 293T , HeLa , HT-29 and HEPG2 cell lysate
<b>WB Recommended dilution:</b>	200-1000

