

## 兔抗 SMC5 多克隆抗体

中文名称：兔抗 SMC5 多克隆抗体

英文名称：Anti-SMC5 rabbit polyclonal antibody

级别：SMC5L1

储存：冷冻（-20℃）

抗原：SMC5

宿主：Rabbit

相关类别：一抗

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

**Background:**

Core component of the SMC5-SMC6 complex, a complex involved in repair of DNA double-strand breaks by homologous recombination. The complex may promote sister chromatid homologous recombination by recruiting the SMC1-SMC3 cohesin complex to double-strand breaks. The complex is required for telomere maintenance via recombination in ALT (alternative lengthening of telomeres) cell lines and mediates sumoylation of shelterin complex (telosome) components which is proposed to lead to shelterin complex disassembly in ALT-associated PML bodies (APBs). Required for recruitment of telomeres to PML nuclear bodies. Required for si

	ster chromatid cohesion during prometaphase and mitotic progression; the function seems to be independent of SMC6.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	SMC5
<b>Immunogen:</b>	Synthetic peptide of human SMC5
<b>Full name:</b>	structural maintenance of chromosomes 5
<b>Synonyms:</b>	SMC5L1
<b>SwissProt:</b>	Q8IY18
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human colorectal cancer
<b>IHC Recommend dilution:</b>	25-100

