

## 兔抗 SMOC1 多克隆抗体

中文名称: 兔抗 SMOC1 多克隆抗体

英文名称: Anti-SMOC1 rabbit polyclonal antibody

别名: SPARC related modular calcium binding 1; OAS

相关类别: 一抗

储存: 冷冻 (-20℃)

抗原: SMOC1

宿主: Rabbit

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a multi-domain secreted protein that may have a critical role in ocular and limb development. Mutations in this gene are associated with microphthalmia and limb anomalies. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	SMOC1
<b>Immunogen:</b>	Fusion protein of human SMOC1
<b>Full name:</b>	SPARC related modular calcium binding 1
<b>Synonyms:</b>	OAS
<b>SwissProt:</b>	Q9H4F8

<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human lung cancer
<b>IHC Recommend dilution:</b>	40-200

