

兔抗 SMUG1 多克隆抗体

中文名称：兔抗 SMUG1 多克隆抗体

英文名称：Anti-SMUG1 rabbit polyclonal antibody

别名：FDG; UNG3; HMUDG

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：SMUG1

反应种属：Human

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	This gene encodes a protein that participates in base excision repair by removing uracil from single- and double-stranded DNA. Many alternatively spliced transcript variants exist for this gene; the full-length nature is known for some but not all of the variants
Applications:	WB
Name of antibody:	SMUG1
Immunogen:	Synthesized peptide derived from internal of human SMUG1.
Full name:	single-strand-selective monofunctional uracil-DNA glycosylase 1
Synonyms :	FDG; UNG3; HMUDG

SwissProt:	Q53HV7
WB Predicted band size:	30 kDa
WB Positive control:	293 cells and COS7 cells lysates
WB Recommended dilution:	500-3000

