

兔抗 SNCB 多克隆抗体

中文名称: 兔抗 SNCB 多克隆抗体

英文名称: Anti-SNCB rabbit polyclonal antibody

相关类别: 一抗

储 存: 冷冻 (-20℃)

宿 主: Rabbit

抗 原: SNCB

反应种属: Human, Mouse, Rat

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a member of a small family of proteins that inhibit phospholipase D2 and may function in neuronal plasticity. The encoded protein is abundant in lesions of patients with Alzheimer disease. A mutation in this gene was found in individuals with dementia with Lewy bodies. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	SNCB
Immunogen:	Synthetic peptide of human SNCB
Full name:	synuclein beta
SwissProt:	Q16143
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer and human tonsil

IHC Recommend dilution:

20-100



