

兔抗 NDUFA12 多克隆抗体

中文名称：兔抗 NDUFA12 多克隆抗体

英文名称：Anti-NDUFA12 rabbit polyclonal antibody

别名：B17.2; DAP13

抗原：NDUFA12

相关类别：一抗

储存：冷冻 (-20℃)

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a protein which is part of mitochondrial complex 1, part of the oxidative phosphorylation system in mitochondria. Complex 1 transfers electrons to ubiquinone from NADH which establishes a proton gradient for the generation of ATP. Mutations in this gene are associated with Leigh syndrome due to mitochondrial complex 1 deficiency. Pseudogenes of this gene are located on chromosomes 5 and 13. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	NDUFA12
Immunogen:	Full length fusion protein
Full name:	NADH dehydrogenase (ubiquinone) 1 alpha subcomplex, 12

Synonyms :	B17.2; DAP13
SwissProt:	Q9UI09
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human lung cancer
IHC Recommend dilution:	25-100

