

兔抗 SNRPD1 多克隆抗体

中文名称: 兔抗 SNRPD1 多克隆抗体

英文名称: Anti-SNRPD1 rabbit polyclonal antibody

别名: small nuclear ribonucleoprotein D1 polypeptide; SMD1; SNRPD; Sm-D1; HsT2456

储存: 冷冻 (-20℃)

抗原: SNRPD1

宿主: Rabbit

反应种属: Human, Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a small nuclear ribonucleoprotein that belongs to the SNRNP core protein family. The protein may act as a charged protein scaffold to promote SNRNP assembly or strengthen SNRNP-SNRNP interactions through nonspecific electrostatic contacts with RNA. Two transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	SNRPD1
Immunogen:	Fusion protein of human SNRPD1
Full name:	small nuclear ribonucleoprotein D1 polypeptide

Synonyms:	SMD1; SNRPD; Sm-D1; HsT2456
SwissProt:	P62314
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and Human cervical cancer
IHC Recommend dilution:	30-150



