

兔抗 SPINK7 多克隆抗体

中文名称: 兔抗 SPINK7 多克隆抗体

英文名称: Anti-SPINK7 rabbit polyclonal antibody

别名: ECG2; ECRG2

抗原: SPINK7

储存: 冷冻 (-20℃)

宿主: Rabbit

相关类别: 一抗

反应种属: Human

标记物: Unconjugate


克隆类型: rabbit polyclonal

技术规格

Background:

SPINK7 (serine peptidase inhibitor, Kazal type 7), also known as Ecg2 (esophagus cancer-related gene 2 protein), is a 76 amino acid secreted protein. Containing one Kazal-like domain, SPINK7 is thought to be a serine protease inhibitor. The gene that encodes SPINK7 maps to human chromosome 5, which contains 181 million base pairs encoding around 1,000 genes. Chromosome 5 is associated with Cockayne syndrome through the ERCC8 gene and familial adenomatous polyposis through the adenomatous polyposis coli (APC) tumor suppressor gene. Treacher Collins syndrome is also chromosome 5 associated and is caused by insertions or

	deletions within the TCOF1 gene. Deletion of the p arm of chromosome 5 leads to Cri du chat syndrome. Deletion of 5q or chromosome 5 altogether is common in therapy-related acute myelogenous leukemias and myelodysplastic syndrome. Probable serine protease inhibitor.
Applications:	ELISA, IHC
Name of antibody:	SPINK7
Immunogen:	Synthetic peptide of human SPINK7
Full name:	serine peptidase inhibitor, Kazal type 7 (putative)
Synonyms :	ECG2; ECRG2
SwissProt:	P58062
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human prostate cancer and Human lung cancer
IHC Recommend dilution:	25-100



mlbio 梅联生物
Good elisakit producers