

## 兔抗 SPIRE1 多克隆抗体

中文名称：兔抗 SPIRE1 多克隆抗体

英文名称：Anti-SPIRE1 rabbit polyclonal antibody

别名：spire type actin nucleation factor 1; Spir-1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：SPIRE1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Spire proteins, such as SPIRE1, are highly conserved between species. They belong to the family of Wiskott-Aldrich homology region-2 (WH2) proteins, which are involved in actin organization (Kerckhoff et al., 2001 [PubMed 11747823]).
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	SPIRE1
<b>Immunogen:</b>	Fusion protein of human SPIRE1
<b>Full name:</b>	spire type actin nucleation factor 1
<b>Synonyms:</b>	Spir-1
<b>SwissProt:</b>	Q08AE8
<b>ELISA Recommended dilution:</b>	5000-10000

<b>IHC positive control:</b>	Human kidney cancer and Human brain
<b>IHC Recommend dilution:</b>	50-300



