

兔抗 SPRTN 多克隆抗体

中文名称: 兔抗 SPRTN 多克隆抗体

英文名称: Anti-SPRTN rabbit polyclonal antibody

别名: DVC1; PRO4323; spartan; C1orf124

储存: 冷冻 (-20℃)

抗原: SPRTN

宿主: Rabbit

相关类别: 一抗

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:

The protein encoded by this gene may play a role in DNA repair during replication of damaged DNA . This protein recruits valosin containing protein (p 97) to stalled DNA replication forks where it may prevent excessive translesional DNA synthesis and limit the number of DNA-damage induced mutations. It may also be involved in replication-related G2/M-checkpoint regulation. Deficiency of a similar protein in mouse causes chromosomal instability and progeroid phenotypes. Mutations in this gene have been associated with Ruijs-Aalfs syndrome (R JALS). Alternatively spliced transcript variants have

	been identified.
Applications:	ELISA, IHC
Name of antibody:	SPRTN
Immunogen:	Fusion protein of human SPRTN
Full name:	SprT-like N-terminal domain
Synonyms:	DVC1; PRO4323; spartan; C1orf124
SwissProt:	Q9H040
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer and human ovarian cancer
IHC Recommend dilution:	50-300



