

兔抗 NID2 多克隆抗体

中文名称: 兔抗 NID2 多克隆抗体

英文名称: Anti-NID2 rabbit polyclonal antibody

别名: nidogen 2; NID-2

储存: 冷冻 (-20℃)

抗原: NID2

宿主: Rabbit

反应种属: Human, Mouse, Rat

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a member of the nidogen family of basement membrane proteins. This protein is a cell-adhesion protein that binds collagens I and IV and laminin and may be involved in maintaining the structure of the basement membrane.
Applications:	ELISA, IHC
Name of antibody:	NID2
Immunogen:	Fusion protein of human NID2
Full name:	nidogen 2
Synonyms:	NID-2
SwissProt:	Q14112
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer

IHC Recommend dilution:	20-100
-------------------------	--------

