

兔抗 NIFK 多克隆抗体

中文名称：兔抗 NIFK 多克隆抗体

英文名称： Anti-NIFK rabbit polyclonal antibody

别名： nucleolar protein interacting with the FHA domain of MKI67; Nopp34; MKI67IP

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： NIFK

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a protein that interacts with the forkhead-associated domain of the Ki-67 antigen. The encoded protein may bind RNA and may play a role in mitosis and cell cycle progression. Multiple pseudogenes exist on chromosomes 5, 10, 12, 15, and 19.
Applications:	ELISA, WB, IHC
Name of antibody:	NIFK
Immunogen:	Fusion protein of human NIFK
Full name:	nucleolar protein interacting with the FHA domain of MKI67
Synonyms:	Nopp34; MKI67IP
SwissProt:	Q9BYG3
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human esophagus cancer and Human thyroid cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	34 kDa
WB Positive control:	A172 cell lysate
WB Recommended dilution:	200-1000



