

兔抗 NLGN1 多克隆抗体

中文名称: 兔抗 NLGN1 多克隆抗体

英文名称: Anti-NLGN1 rabbit polyclonal antibody

别名: NL1

储存: 冷冻 (-20℃)

抗原: NLGN1

宿主: Rabbit

相关类别: 一抗

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a member of a family of neuronal cell surface proteins. Members of this family may act as splice site-specific ligands for beta-neurexins and may be involved in the formation and remodeling of central nervous system synapses.
Applications:	ELISA, IHC
Name of antibody:	NLGN1
Immunogen:	Synthetic peptide of human NLGN1
Full name:	neuroligin 1
Synonyms:	NL1
SwissProt:	Q8N2Q7
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer

IHC Recommend dilution: 20-100

