

兔抗 STEAP2 多克隆抗体

- 中文名称: 兔抗 STEAP2 多克隆抗体
- 英文名称: Anti-STEAP2 rabbit polyclonal antibody
- 别名: STMP; IPCA1; PUMPCn; STAMP1; PCANAP1
- 相关类别: 一抗
- 储存: 冷冻 (-20℃)
- 宿主: Rabbit
- 抗原: STEAP2
- 反应种属: Human, Mouse
- 标记物: Unconjugate
- 克隆类型: Unconjugate

技术规格

Background:	This gene is a member of the STEAP family and encodes a multi-pass membrane protein that localizes to the Golgi complex, the plasma membrane, and the vesicular tubular structures in the cytosol. A highly similar protein in mouse has both ferrireductase and cupric reductase activity, and stimulates the cellular uptake of both iron and copper in vitro. Increased transcriptional expression of the human gene is associated with prostate cancer progression. Alternate transcriptional splice variants, encoding different isoforms, have been characterized.
Applications:	ELISA, WB

Name of antibody:	STEAP2
Immunogen:	Synthetic peptide of human STEAP2
Full name:	STEAP family member 2, metalloredutase
Synonyms :	STMP; IPCA1; PUMPCn; STAMP1; PCANAP1
SwissProt:	Q8NFT2
ELISA Recommended dilution:	5000-10000
WB Predicted band size:	56 kDa
WB Positive control:	Human prostate tissue
WB Recommended dilution:	500-2000

