

兔抗 STK17A 多克隆抗体

中文名称: 兔抗 STK17A 多克隆抗体

英文名称: Anti-STK17A rabbit polyclonal antibody

别名: DRAK1

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: STK17A

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene is a member of the DAP kinase-related apoptosis-inducing protein kinase family and encodes an autophosphorylated nuclear protein with a protein kinase domain. The protein has apoptosis-inducing activity.
Applications:	ELISA, IHC
Name of antibody:	STK17A
Immunogen:	Synthetic peptide of human STK17A
Full name:	serine/threonine kinase 17a
Synonyms:	DRAK1
SwissProt:	Q9UEE5
IHC positive control:	Human gastric cancer and human brain
IHC Recommend dilution:	25-100

