

## 兔抗 STK32B 多克隆抗体

中文名称: 兔抗 STK32B 多克隆抗体

英文名称: Anti-STK32B rabbit polyclonal antibody

别名: STK32; STKG6; YANK2; HSA250839

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: STK32B

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	STK32B is a 414 amino acid protein that contains one protein kinase domain and belongs to the Ser/Thr protein kinase family. Using magnesium as a cofactor, YANK2 functions to catalyze the ATP-dependent phosphorylation of target proteins and may play a role in various signaling events throughout the cell. Multiple isoforms of YANK2 exist due to alternatively splicing events.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	STK32B
<b>Immunogen:</b>	Synthetic peptide of human STK32B
<b>Full name:</b>	serine/threonine kinase 32B

<b>Synonyms :</b>	STK32; STKG6; YANK2; HSA250839
<b>SwissProt:</b>	Q9NY57
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human cervical cancer and Human prostate cancer
<b>IHC Recommend dilution:</b>	25-100

