

兔抗 NR2F1 多克隆抗体

中文名称：兔抗 NR2F1 多克隆抗体

英文名称：Anti-NR2F1 rabbit polyclonal antibody

别名：nuclear receptor subfamily 2 group F member 1 EAR3; BBOAS; EAR-3; NR2F2; SVP44; BBSOAS; ERBAL3; TFCOUP1; COUP-TFI; TCFCOUP1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：NR2F1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a nuclear hormone receptor and transcriptional regulator. The encoded protein acts as a homodimer and binds to 5'-AGGTCA-3' repeats. Defects in this gene are a cause of Bosch-Boonstra optic atrophy syndrome (BBOAS).
Applications:	ELISA, WB, IHC
Name of antibody:	NR2F1
Immunogen:	Synthetic peptide of human NR2F1
Full name:	nuclear receptor subfamily 2 group F member 1
Synonyms:	EAR3; BBOAS; EAR-3; NR2F2; SVP44; BBSOAS; ERBAL3; TFCOUP1; COUP-TFI; TCFCOUP1

SwissProt:	P10589
ELISA Recommended dilution :	5000-10000
IHC positive control:	Human ovarian cancer
IHC Recommend dilution:	20-100
WB Predicted band size:	46 kDa
WB Positive control:	293T and NIH/3T3 cell lysates
WB Recommended dilution:	200-1000



