

## 兔抗 STRAP 多克隆抗体

中文名称: 兔抗 STRAP 多克隆抗体

英文名称: Anti-STRAP rabbit polyclonal antibody

别名: MAWD; PT-WD; UNRIP

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: STRAP

相关类别: 一抗

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	Serine-threonine kinase receptor-associated protein is an enzyme that in humans is encoded by the STRAP gene. The SMN complex plays a catalyst role in the assembly of small nuclear ribonucleoproteins (snRNPs), the building blocks of the spliceosome. Thereby, plays an important role in the splicing of cellular pre-mRNAs. Most spliceosomal snRNPs contain a common set of Sm proteins SNRPB, SNRPD1, SNRPD2, SNRPD3, SNRPE, SNRPF and SNRPG that assemble in a heptameric protein ring on the Sm site of the small nuclear RNA to form the core snRNP.
<b>Applications:</b>	ELISA, IHC

<b>Name of antibody:</b>	STRAP
<b>Immunogen:</b>	Full length fusion protein
<b>Full name:</b>	serine/threonine kinase receptor associated protein
<b>Synonyms :</b>	MAWD; PT-WD; UNRIP
<b>SwissProt:</b>	Q9Y3F4
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human lung cancer
<b>IHC Recommend dilution:</b>	20-100

