

## 兔抗 SUN5 多克隆抗体

中文名称：兔抗 SUN5 多克隆抗体

英文名称：Anti-SUN5 rabbit polyclonal antibody

别名：SPAG4L; TSARG4; dJ726C3.1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：SUN5

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

**Background:**

TSARG4 (testis and spermatogenesis-related gene 4 protein), also known as SPAG4L (sperm-associated antigen 4-like protein) or SRG4 (spermatogenesis-related gene 4) in mice, is a ubiquitously expressed protein in humans but, in mice, it is specifically expressed in the testis. TSARG4 belongs to the spindle pole body family of proteins and contains one Unc84/SUN domain and a potential transmembrane region. The human TSARG4 protein shares 74% amino acid identity with its mouse counterpart. In mice, TSARG4 is found in round spermatids and spermatocytes and may participate in mitosis. It has been hypothesized th

	at TSARG4 localizes to the nuclear envelope.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	SUN5
<b>Immunogen:</b>	Synthetic peptide of human SUN5
<b>Full name:</b>	Sad1 and UNC84 domain containing 5
<b>Synonyms:</b>	SPAG4L; TSARG4; dJ726C3.1
<b>SwissProt:</b>	Q8TC36
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human lung cancer and human brain
<b>IHC Recommend dilution:</b>	20-100



