

兔抗 NSD2 多克隆抗体

中文名称：兔抗 NSD2 多克隆抗体

英文名称：Anti-NSD2 rabbit polyclonal antibody

别名：nuclear receptor binding SET domain protein 2; WHS

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：NSD2

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a protein that contains four domains present in other developmental proteins: a PWWP domain, an HMG box, a SET domain, and a PHD-type zinc finger. It is expressed ubiquitously in early development. Wolf-Hirschhorn syndrome (WHS) is a malformation syndrome associated with a hemizygous deletion of the distal short arm of chromosome 4. This gene maps to the 165 kb WHS critical region and has also been involved in the chromosomal translocation t(4;14)(p16.3;q32.3) in multiple myelomas. Alternative splicing of this gene results in multiple transcript variants encoding different isoforms. Some transcript variants are nonsense-mediated mRNA (NMD) decay candidates, hence not represented as reference sequences. [provided by RefSeq, Jul 2008]
Applications:	ELISA, IHC

Name of antibody:	NSD2
Immunogen:	Synthetic peptide of human NSD2
Full name:	nuclear receptor binding SET domain protein 2
Synonyms:	WHS; TRX5; KMT3F; KMT3G; MMSET; WHSC1; REIIBP
SwissProt:	O96028
IHC positive control:	Human colorectal cancer
IHC Recommend dilution:	25-100

