

## 兔抗 NTN1 多克隆抗体

中文名称：兔抗 NTN1 多克隆抗体

英文名称：Anti-NTN1 rabbit polyclonal antibody

别名：NTN1L

抗原：NTN1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

<b>Background:</b>	Netrin is included in a family of laminin-related secreted proteins. The function of this gene has not yet been defined; however, netrin is thought to be involved in axon guidance and cell migration during development. Mutations and loss of expression of netrin suggest that variation in netrin may be involved in cancer development.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	NTN1
<b>Immunogen:</b>	Synthetic peptide of human NTN1
<b>Full name:</b>	netrin 1
<b>Synonyms :</b>	NTN1L
<b>SwissProt:</b>	O95631
<b>ELISA Recommended dilution:</b>	5000-10000

<b>IHC positive control:</b>	Human tonsil
<b>IHC Recommend dilution:</b>	25-100

