

## 兔抗 TAC1 多克隆抗体

中文名称：兔抗 TAC1 多克隆抗体

英文名称：Anti-TAC1 rabbit polyclonal antibody

别名：NK2, NPK, NKNA, TAC2, Hs.2563

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：TAC1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes four products of the tachykinin peptide hormone family, substance P and neurokinin A, as well as the related peptides, neuropeptide K and neuropeptide gamma. These hormones are thought to function as neurotransmitters which interact with nerve receptors and smooth muscle cells. They are known to induce behavioral responses and function as vasodilators and secretagogues. Multiple transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	TAC1
<b>Immunogen:</b>	Synthetic peptide of human TAC1
<b>Full name:</b>	Tachykinin, precursor 1

<b>Synonyms :</b>	NK2, NPK, NKNA, TAC2, Hs.2563
<b>SwissProt:</b>	P20366
<b>ELISA Recommended dilution:</b>	1000-5000
<b>IHC positive control:</b>	Human esophagus cancer
<b>IHC Recommend dilution:</b>	25-100

